

wee SECURE™

Instructions

Intended Use:

The wee SECURE is a reusable multi-patient device intended for the temporary immobilization of infants during circumcision and other procedures, such as placement of intravenous infusions, x-rays, and minor procedures

Specifications

The wee SECURE fits infants weighing 1800 - 9072 grams (3.9 - 20 pounds).

Using the wee SECURE



CAUTION


- Read instructions completely before using this product.
- This product should be used only under the direct supervision of a licensed medical practitioner. "Rx Only."
- While immobilized, patient should be under continuous supervision. **Never leave patient unattended.**
- Board to be used in standard incubator or bassinett.

Required Items

- wee SECURE
- wee STRAPS (one for each arm and leg)

To Install the Straps:

- 1) Place board upside down without infant.
- 2) Thread the arm and leg straps through the desired slots in the direction of intended use.



Product Material


- Board: Dense, nonabsorbant plastic - contains no metal, MR safe
- Straps: Polyurethane Foam/Polyester Fleece Tricot
- wee SECURE and wee STRAPS are latex-free

Product Size

- Board: 14 3/4" (37.46 cm) wide x 25" (63.50 cm) long x 2" (5.08 cm) deep, 1.85 lbs (839 g)
- Straps: 11.4 " (28.95 cm) each, .3 oz (85 g)

To Immobilize the Infant:

- 1) A barrier should be used between infant and board.
- 2) Place the infant in the board. Position the straps as near as possible to the elbows/arms and knees.
- 3) Fasten the straps by gently pressing one end of the strap on the infant while bringing the opposite end of the strap across the top and snugly fasten the hooks and loops of the strap together.



To Position for Procedure/Circumcision:

The infant should be carefully placed into wee SECURE so that the body is inside the depressions of the board.

- * Position and fasten the wee STRAPS
- * Use discretion for strap placement based on size of patient and type of procedure.
- * Towels and/or padding should be used for additional comfort of the infant.

To Clean the wee SECURE Board:

- 1) Clean board with soap and water.
- 2) Disinfect using medical-grade disinfectant.
- 3) Rinse board and dry completely.



CAUTION

Do not heat sterilize the board. Do not clean with agents containing ketones, ether, or other plastic solvents.

Storing the wee SECURE

Hang on a wall using the hole on side of the board near head or nest several wee SECUREs together on a shelf.

Required Items

- Disinfectant liquid or wipes
- Dry towel to remove excess disinfectant

Ordering Information

wee SECURE.....2017-1
wee STRAPS.....2017-4
(set of four arm and/or leg straps)

Warranty

Unconditional one-year warranty against all manufacturing defects.

wee STRAPS available at: