<span><span><span></span></span></span> <div>©2019 wee Medical</div>
<p>קילפא צװישן די גלידער "wee BELL" רינגל די דיגריטן</p>
<p>ukubwa sita.</p> <p>Kifaa cha tohara cha Weusi kinakujia kwa</p>
<p>El pequeño dispositivo de circuncisión Bell viene en seis tamaños.</p>
<p>Das kleine BELL-Beschneidungsgerät ist in sechs Größen erhältlich.</p>
<p>Le dispositif de circoncision wee BELL est disponible en six tailles.</p>

In six sizes.

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
**Using the wee BELL Circumcision Device**

**Caution:** The glans must be free in its entirety under the device.

Pressure is caused by 1) pulling the foreskin too far toward in an attempt to remove maximum tissue, or 2) using a device that is too small.

Using a device that is too large and/or pulling the foreskin too far forward may cause the BELL to place pressure on the glans and meatus. Edema may result and contribute to damage to the skin and underlying tissue. It would then be necessary to remove the device by cutting the ligature and pulling off the BELL.

The ligature must be tied securely in the groove with a square knot.

Parents should be instructed to notify the care provider if the wee BELL has not fallen off within eight days, or immediately if it cannot be exposed or inspected or if the penis becomes edematous. Information leaflets provided should be given to parents to supplement the care provider’s instructions.

**How to Use the wee BELL Circumcision Device**

**1.** Apply two hemostats to foreskin, and gently separate the adhesions using either the pointed or spatula end of the probe. Use the hemostats to spread the foreskin and another to crush the foreskin where the dorsal slit is to be made. The foreskin should be crushed to a length that is approximately the same as the width of the glans. This can be visualized while the foreskin is being pulled and dissected. Place ligature with the beginning of a square knot or surgeon’s knot and leave loose at the base of the penis.

**2.** With tissue scissors, make the dorsal slit along the crushed line. The slit should be no longer than necessary to permit the BELL to be worked into place, and it should not exceed the length of the previously crushed area. Gently retract the foreskin and free any remaining adhesions, completely exposing sulcus.

**3.** Place BELL of appropriate size (see below) over the glans so it avoids undue pressure on frenulum area. Too small a BELL may injure the glans. Too large a BELL may cause BELL to shift, possibly causing edema and damage to skin and underlying tissue. Positioning the wee BELL so that the groove on the handle is aligned with the urethra will help prevent BELL pressure on the urethra in case the BELL shifts while in place. Apex of dorsal slit should be distal to BELL’s groove. Pull foreskin enough for the desired outcome but always avoid positioning apex of dorsal slit proximal to groove. If excessively pulled, it may later draw the BELL too tightly against the glans onto the shaft. There should be approximately 1.0 cm between the sulcus and the groove of the BELL where the ligature will be applied.

Additional suggestions:

- To hold BELL while tying ligature, a third hemostat may be used to clamp foreskin to BELL handle or a safety pin can be passed through the hole in the BELL handle.
- If slit is too long, cross the two hemostats over each other to hold the BELL in place.

**4.** After positioning ligature around and into the BELL’s groove, draw ligature very tightly so it compresses foreskin into groove; tie with a square knot or surgeon’s knot. Cut off excess ligature leaving approximately 5mm length. You should be able to see an unobstruct- ed urethral meatus if the groove on the handle is in alignment with the urethra.

**5.** Trim off foreskin using outer ridge of BELL as cutting guide. Maintain adequate traction on foreskin and trim as much tissue as possible to reduce the amount of remaining necrotic tissue. If ligature is not tightly secured, bleeding may occur and BELL may fall off prematurely, or the line of excision will be incomplete, requiring subsequent removal of tissue tags.

**6.** Break off BELL handle and discard. Leaving the BELL and ligature in place, you should be able to see an unobstructed urethral meatus. No dressing is necessary. The rim of tissue under and distal to the ligature will become necrotic and separate with the BELL in 2 to 8 days, leaving a clean, healed line of excision. If the care provider desires, the BELL can be removed before the infant is sent home by cutting the ligature at the knot. This may be done after 36 to 48 hours following circumcision, based on care provider’s evaluation.

**Proper Size Selection**

Proper size selection is most important since the wee BELL device may remain in place for as long as eight days. After the foreskin has been retracted (Step 2), the proper wee BELL size can easily be determined.

Too Small: wee BELL device too small may cause the ring to exert undue pressure on the sides of the glans. Use of a BELL this size may cause local necrotic changes and periurethral edema to the point of causing urinary retention. A too-small BELL may also result in insufficient circumcision.

Too Large: wee BELL device too large may shift and torque and cause pressure on penis shaft or urethra. Use of a too-large BELL may cause over-stretching of the foreskin and increase the possibility of edema. Too much foreskin may be removed using the BELL this size.

**Warning:** To reduce the potential for infection do not re-use.

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